

Efficient

Effective

Economical



8800 shown with nebulizer kit attached



8800 shown


Rescuer[®] II Compact CPAP System

- **Lowest Oxygen Consumption will help your tank last twice as long as other CPAP system, which will reduce your on-site call times.**
- **Only 8 LPM to achieve 10 cmH₂O**
- **Pre-attached filter to reduce spread of possible contaminants**
- **Single flowmeter required to provide CPAP and nebulized medication**
- **Rapid deployment using high quality CPAP mask with adjustable forehead pads**
- **Available with nebulizer kit with 24" oxygen tubing and T-piece**



The Rescuer Emergency CPAP System developed by BLS Systems is a state of the art non-invasive ventilator support system. Controlled airway pressure provides rapid relief for maximum patient benefit with minimal oxygen consumption

Performance - 8Lpm delivers 10cmH2O



CPAP Setting	5cmH2O	5cmH2O	6cmH2O	7.5cmH2O	10cmH2O	15cmH2O
Set Oxygen Flow	5Lpm	6Lpm	7Lpm	8Lpm	9Lpm	10Lpm
Oxygen Percent	53%	52%	52%	53%	57%	60%

Product Information

Product Number	Description
8800	RESCUER® II Compact CPAP System with adjustable Large CPAP Mask
8805	RESCUER® II Compact CPAP System with adjustable Medium CPAP Mask
8810	RESCUER® II Compact CPAP System with adjustable Small CPAP Mask
8850	RESCUER®II Compact CPAP System with #6 Cushion Mask
8855	RESCUER® II Compact CPAP System with #4 Cushion Mask
8900	Nebulizer Kit for RESCUER®CPAP Systems (nebulizer, T-piece, 24" tubing)
8800-NK	8800 System with Nebulizer Kit
8805-NK	8805 System with Nebulizer Kit
8810-NK	8810 System with Nebulizer Kit

BLS Systems Limited
1124 South Service Road West
Oakville, Ontario
Canada L6L 5T7

Telephone: (905)-339-1069
Fax: (905)-339-1060
E-mail: sales@blssystemsltd.com
Website: http://www.blssystemsltd.com