



Rescuer[®] II Compact CPAP System

Efficient



8800 shown

Effective




8850 shown

- Low cost, completely disposable and simple to use
- Highly efficient for longer cylinder life – runs between 4 to 8Lpm
- Lowest oxygen consumption of any disposable CPAP device
- Maximum protection - provides filtration on inhalation AND exhalation
- Allows medication delivery
- Simple operation
- Medication Port
- Filters (3M HEPA)
- High quality CPAP Mask with adjustable forehead pads (8800, 8805)
- Integral manometer
- Critical filter placement does not separate patient from CPAP generator
- Deliver CPAP with nebulized medications using a **single** oxygen source.
- Patent Pending

Rescuer[®] II Compact CPAP System

The Rescuer II Compact CPAP System developed by BLS Systems is a state of the art non-invasive ventilator support system. Controlled airway pressure provides rapid relief for maximum patient benefit with minimal oxygen consumption

Performance - Achieve 15cmH₂O at 8 Lpm



Product Information

Product Number	Description
8800	Rescuer [®] II Compact CPAP System with adjustable Large CPAP Mask
8805	Rescuer [®] II Compact CPAP System with adjustable Medium CPAP Mask
8850	Rescuer [®] II Compact CPAP System with #6 Cushion Mask
8855	Rescuer [®] II Compact CPAP System with #4 Cushion Mask

BLS Systems Limited
1124 South Service Road West
Oakville, Ontario
Canada L6L 5T7

Telephone: (905)-339-1069
Fax: (905)-339-1060
E-mail: sales@blssystemsltd.com
Website: <http://www.blssystemsltd.com>